



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ANKRD20A3 Antibody [#C14494]

Catalog Number: C14494

Concentration: 1mg/ml

Swiss-Prot No.: Q5VUR7

Other Names: Ankyrin repeat domain-containing protein 20A3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

ANKRD20A3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity: ANKRD20A3 antibody detects endogenous levels of total ANKRD20A3 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IHC: 1:50~1:100 ELISA: 1:10000

References:

Humphray S.J., Nature 429:369-374(2004).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human brain tissue using ANKRD20A3 antibody.