



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ANKRD26 Antibody [#C14497]

Catalog Number: C14497

Concentration: 1mg/ml

Swiss-Prot No.: Q9UPS8

Other Names: ankyrin repeat domain 26; ANR26

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ANKRD26.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ANKRD26 antibody detects endogenous levels of total ANKRD26 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)


Applications: IHC: 1:50~1:100 ELISA: 1:5000

References:

Kikuno R., DNA Res. 6:197-205(1999).

Deloukas P., Nature 429:375-381(2004).

Ota T., Nat. Genet. 36:40-45(2004).



Peptide

-

+

Immunohistochemistry analysis of paraffin-embedded human brain tissue using ANKRD26 antibody.