



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ANXA10 Antibody [#C14499]

Catalog Number: C14499

Concentration: 1mg/ml

Swiss-Prot No.: Q9UJ72

Other Names: annexin 14; annexin A10; annexin-14; ANX10; ANX14

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human ANXA10.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: ANXA10 antibody detects endogenous levels of total ANXA10 protein.

Reactivity: Human (Identities = 100%, Positives = 100%);
Mouse (Identities = 92%, Positives = 100%)

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Morgan R.O., Genomics 60:40-49(1999).

Paul H., Submitted (OCT-1999) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).

