



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

APOL1 Antibody [#C14539]

Catalog Number: C14539

Concentration: 1mg/ml

Swiss-Prot No.: O14791

Other Names: Apo-L; APOL; ApoL-I; Apolipoprotein L;
Apolipoprotein L-I; Apolipoprotein-L1 precursor

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human APOL1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: APOL1 antibody detects endogenous levels of total APOL1 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Duchateau P.N., J. Biol. Chem. 272:25576-25582(1997).

Duchateau P.N., J. Lipid Res. 42:620-630(2001).

Page N.M., Genomics 74:71-78(2001).

