



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

AQP12 Antibody [#C14549]

Catalog Number: C14549

Concentration: 1mg/ml

Swiss-Prot No.: Q8IXF9

Other Names: Aquaporin-12A; AQP-12; AQP12A; AQP12, AQPX2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human AQP12.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: AQP12 antibody detects endogenous levels of total AQP12 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Itoh T., Biochem. Biophys. Res. Commun. 330:832-838(2005).

Ishibashi K., Publishers, New York. (2000).

