



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ARG1 Antibody [#C14559]

Catalog Number: C14559

Concentration: 1mg/ml

Swiss-Prot No.: P05089

Other Names: AI; AI256583; Arg-1; ARG11; arginase 1, liver; arginase I; PGIF

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ARG1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ARG1 antibody detects endogenous levels of total ARG1 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)


Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Haraguchi Y., Proc. Natl. Acad. Sci. U.S.A. 84:412-415(1987).

Takiguchi M., Nucleic Acids Res. 16:8789-8802(1988).

Mungall A.J., Nature 425:805-811(2003).



For Research Use Only