



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

ARMC6 Antibody [#C14562]

Catalog Number: C14562

Concentration: 1mg/ml

Swiss-Prot No.: Q6NXE6

Other Names: armadillo repeat containing 6

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human ARMC6.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: ARMC6 antibody detects endogenous levels of total ARMC6 protein.

Reactivity: Human, Mouse, Rat


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000

References:


Grimwood J., Nature 428:529-535(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Matsuoka S., Science 316:1160-1166(2007).



Western blot analysis of extracts from HeLa cells and COLO cells, using ARMC6 antibody.



Immunofluorescence analysis of A549 cells, using ARMC6 antibody.