



Web: www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR51T1 Antibody [#G455]

Catalog Number: G455

Concentration: 1mg/ml

Swiss-Prot No.: Q8NGJ9

Other Names: Olfactory receptor 51T1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human OR51T1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity: OR51T1 antibody detects endogenous levels of total OR51T1 protein.

Reactivity: Human


Applications: WB: 1:500~1:1000 IF: 1:100~1:500
ELISA: 1:10000

References:

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.



Western blot analysis of extracts from HT-29 cells and A549 cells, using OR51T1 antibody.



Immunofluorescence analysis of A549 cells, using OR51T1 antibody.