



Web: www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR52A4 Antibody [#G456]

Catalog Number: G456

Concentration: 1mg/ml

Swiss-Prot No.: A6NMU1

Other Names: Putative olfactory receptor 52A4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

OR52A4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR52A4 antibody detects endogenous levels of total OR52A4 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Taylor T.D., Nature 440:497-500(2006).

