



Web: www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR52E6 Antibody [#G462]

Catalog Number: G462

Concentration: 1mg/ml

Swiss-Prot No.: Q96RD3

Other Names: olfactory receptor 52E6; Q96RD3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human OR52E6.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR52E6 antibody detects endogenous levels of total OR52E6 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

Fuchs T., Genomics 80:295-302(2002).

