



Web: www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR52N4 Antibody [#G465]

Catalog Number: G465

Concentration: 1mg/ml

Swiss-Prot No.: Q8NGI2

Other Names: Olfactory receptor 52N4; Olfactory receptor OR11-64; OR52N4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human OR52N4.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR52N4 antibody detects endogenous levels of total OR52N4 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

Malnic B., Proc. Natl. Acad. Sci. U.S.A. 101:2584-2589(2004).

Menashe I., Nat. Genet. 34:143-144(2003).

