



Web: www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR52W1 Antibody [#G466]

Catalog Number: G466

Concentration: 1mg/ml

Swiss-Prot No.: Q6IF63

Other Names: O52W1; olfactory receptor 52W1; olfactory receptor OR11-71; olfactory receptor, family 52, subfamily W, member 1; OR11-71; OR52W1; OR52W1P

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human OR52W1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR52W1 antibody detects endogenous levels of total OR52W1 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

Taylor T.D., Nature 440:497-500(2006).

Malnic B., Proc. Natl. Acad. Sci. U.S.A. 101:2584-2589(2004).

