



Web: www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR56B4 Antibody [#G470]

Catalog Number: G470

Concentration: 1mg/ml

Swiss-Prot No.: Q8NH76

Other Names: O56B4; olfactory receptor 56B4; olfactory receptor OR11-67; olfactory receptor, family 56, subfamily B, 4; OR11-67

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

OR56B4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR56B4 antibody detects endogenous levels of total OR56B4 protein.


Reactivity: Human

Applications: WB: 1:500~1:3000 IF: 1:100~1:500


ELISA: 1:1000

References:

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.



Western blot analysis of extracts from COLO cells, using OR56B4 antibody.



Immunofluorescence analysis of MCF-7 cells, using OR56B4 antibody.