



Web: www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR5AP2 Antibody [#G473]

Catalog Number: G473

Concentration: 1mg/ml

Swiss-Prot No.: Q8NGF4

Other Names: Olfactory receptor 5AP2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human OR5AP2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR5AP2 antibody detects endogenous levels of total OR5AP2 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:1000

References:

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

