



Web: www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR6C68 Antibody [#G480]

Catalog Number: G480

Concentration: 1mg/ml

Swiss-Prot No.: A6NDL8

Other Names: O6C68; olfactory receptor 6C68; olfactory receptor, family 6, subfamily C, member 68

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human OR6C68.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR6C68 antibody detects endogenous levels of total OR6C68 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Scherer S.E., Nature 440:346-351(2006).

Mural R.J., Submitted (JUL-2005) to the EMBL/GenBank/DDBJ databases.

