



Web: www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OPRK1 Antibody [#G491]

Catalog Number: G491

Concentration: 1mg/ml

Swiss-Prot No.: P41145

Other Names: kappa-type opioid receptor; MSL-1; OPRK; OPRK1; ROR-D

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human OPRK1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OPRK1 antibody detects endogenous levels of total OPRK1 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Mansson E., Biochem. Biophys. Res. Commun. 202:1431-1437(1994).

Simonin F., Proc. Natl. Acad. Sci. U.S.A. 92:7006-7010(1995).

Zhu J., Life Sci. 56:PL201-PL207(1995).

