



Web: www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

RHO Antibody [#G495]

Catalog Number: G495

Concentration: 1mg/ml

Swiss-Prot No.: P08100

Other Names: CSNBAD1; MGC138309; MGC138311; OPN2; OPSD; opsin 2, rod pigment; Opsin-2; retinitis pigmentosa 4, autosomal dominant; RHO; Rhodopsin; rhodopsin (opsin 2, rod pigment) (retinitis pigmentosa 4, autosomal dominant); RP4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from Internal of human RHO.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: RHO antibody detects endogenous levels of total RHO protein.

Reactivity: Human; Mouse; Rat


Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Nathans J., Proc. Natl. Acad. Sci. U.S.A. 81:4851-4855(1984).

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

Bechtel S., BMC Genomics 8:399-399(2007).



For Research Use Only