



Web: www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR10A6 Antibody [#G497]

Catalog Number: G497

Concentration: 1mg/ml

Swiss-Prot No.: Q8NH74

Other Names: O10A6; olfactory receptor 10A6; olfactory receptor family 10 subfamily A member 6; olfactory receptor OR11-96; OR11-96

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human OR10A6.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR10A6 antibody detects endogenous levels of total OR10A6 protein.

Reactivity: Human

Applications: WB: 1:500~1:1000 ELISA: 1:1000

References:

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

Malnic B., Proc. Natl. Acad. Sci. U.S.A. 101:2584-2589(2004).

