



Web: [www.dijibio.com](http://www.dijibio.com)  
Order: [order@dijibio.com](mailto:order@dijibio.com)  
Tel: +86 519 8805 0026

## OR10D4 Antibody [#G498]

**Catalog Number:** G498

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q8NGN7

**Other Names:** Putative olfactory receptor 10D4;  
OR10D4P; OR10D4

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human OR10D4.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


**Specificity:** OR10D4 antibody detects endogenous levels of total OR10D4 protein.

**Reactivity:** Human


**Applications:** WB: 1:500~1:3000 IF: 1:100~1:500  
ELISA: 1:10000

**References:**

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.



Western blot analysis of extracts from COLO cells, using OR10D4 antibody.



Immunofluorescence analysis of LOVO cells, using OR10D4 antibody.