



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR10J1 Antibody [#G502]

Catalog Number: G502

Concentration: 1mg/ml

Swiss-Prot No.: P30954

Other Names: HGMP07J; HSHGMP07J; MGC138495; MGC138497; O10J1; olfactory receptor 10J1; olfactory receptor OR1-26; olfactory receptor, family 10, subfamily J, member 1; olfactory receptor-like protein HGMP07J

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

OR10J1.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR10J1 antibody detects endogenous levels of total OR10J1 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Parmentier M., Nature 355:453-455(1992).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Fuchs T., Genomics 80:295-302(2002).

