



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR1L6 Antibody [#G529]

Catalog Number: G529

Concentration: 1mg/ml

Swiss-Prot No.: Q8NGR2

Other Names: HG16; olfactory receptor 1L6; olfactory receptor, family 1, subfamily L, member 6; OR1L7; OR9-30

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human OR1L6.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity: OR1L6 antibody detects endogenous levels of total OR1L6 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:20000

References:

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.
Fuchs T., Genomics 80:295-302(2002).



Western blot analysis of extracts from HUVEC cells and HeLa cells, using OR1L6 antibody.



Immunofluorescence analysis of HUVEC cells, using OR1L6 antibody.