



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR2H2 Antibody [#G551]

Catalog Number: G551

Concentration: 1mg/ml

Swiss-Prot No.: O95918

Other Names: Olfactory receptor 2H2; Olfactory receptor-like protein FAT11; Hs6M1-12; OR2H2; FAT11, OLFR2, OR2H3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

OR2H2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR2H2 antibody detects endogenous levels of total OR2H2 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Fan W., Genomics 27:119-123(1995).

Eklund A.C., Hum. Immunol. 61:711-717(2000).

Mungall A.J., Nature 425:805-811(2003).

