



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR2J3 Antibody [#G554]

Catalog Number: G554

Concentration: 1mg/ml

Swiss-Prot No.: O76001

Other Names: Olfactory receptor 2J3; Olfactory receptor OR6-16; Olfactory receptor 6-6; OR6-6; Hs6M1-3; OR2J3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human OR2J3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR2J3 antibody detects endogenous levels of total OR2J3 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:20000

References:

Mungall A.J., Nature 425:805-811(2003).

Fuchs T., Genomics 80:295-302(2002).

Malnic B., Proc. Natl. Acad. Sci. U.S.A. 101:2584-2589(2004).

