



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR2M2 Antibody [#G560]

Catalog Number: G560

Concentration: 1mg/ml

Swiss-Prot No.: Q96R28

Other Names: olfactory receptor 2M2; olfactory receptor, family 2, subfamily M, member 2; OR2M2Q; OST423

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human OR2M2.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR2M2 antibody detects endogenous levels of total OR2M2 protein.


Reactivity: Human

Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000


References:

Gregory S.G., Nature 441:315-321(2006).

Fuchs T., Genomics 80:295-302(2002).



Western blot analysis of extracts from LOVO cells and A549 cells, using OR2M2 antibody.



Immunofluorescence analysis of MCF-7 cells, using OR2M2 antibody.