



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR2Y1 Antibody [#G574]

Catalog Number: G574

Concentration: 1mg/ml

Swiss-Prot No.: Q8NGV0

Other Names: Olfactory receptor 2Y1; Olfactory receptor OR5-2

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human OR2Y1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR2Y1 antibody detects endogenous levels of total OR2Y1 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000

References:


Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

Fuchs T., Genomics 80:295-302(2002).

Malnic B., Proc. Natl. Acad. Sci. U.S.A. 101:2584-2589(2004).



Western blot analysis of extracts from HT-29 cells and RAW264.7 cells, using OR2Y1 antibody.



Immunofluorescence analysis of A549 cells, using OR2Y1 antibody.