



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** +86 519 8805 0026

## OR3A2 Antibody [#G577]

**Catalog Number:** G577

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P47893

**Other Names:** olfactory receptor OR17-14; olfactory receptor, family 3, subfamily A, member 2; OLFRA04; OR17-14; OR17-228; OR228

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

OR3A2.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** OR3A2 antibody detects endogenous levels of total OR3A2 protein.


**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 IF: 1:100~1:500  
ELISA: 1:10000


**References:**

Glusman G., Genomics 63:227-245(2000).

Ben-Arie N., Hum. Mol. Genet. 3:229-235(1994).



Western blot analysis of extracts from Jurkat cells, using OR3A2 antibody.



Immunofluorescence analysis of MCF-7 cells, using OR3A2 antibody.