



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR4L1 Antibody [#G608]

Catalog Number: G608

Concentration: 1mg/ml

Swiss-Prot No.: Q8NH43

Other Names: olfactory receptor 4L1; OR14-28; OR4L2P

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from C-terminal of human OR4L1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR4L1 antibody detects endogenous levels of total OR4L1 protein.


Reactivity: Human

Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:10000


References:

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

Malnic B., Proc. Natl. Acad. Sci. U.S.A. 101:2584-2589(2004).



Western blot analysis of extracts from K562 cells and MCF-7 cells, using OR4L1 antibody.



Immunofluorescence analysis of COS-7 cells, using OR4L1 antibody.