



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR5AR1 Antibody [#G625]

Catalog Number: G625

Concentration: 1mg/ml

Swiss-Prot No.: Q8NGP9

Other Names: Olfactory receptor 5AR1; Olfactory receptor OR11-209

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human OR5AR1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR5AR1 antibody detects endogenous levels of total OR5AR1 protein.

Reactivity: Human


Applications: WB: 1:500~1:1000 IF: 1:100~1:500
ELISA: 1:10000

References:


Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

Malnic B., Proc. Natl. Acad. Sci. U.S.A. 101:2584-2589(2004).

Menashe I., PLoS Biol. 5:2462-2468(2007).



Western blot analysis of extracts from Jurkat cells, COLO cells and MCF-7 cells, using OR5AR1 antibody.



Immunofluorescence analysis of A549 cells, using OR5AR1 antibody.