



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** +86 519 8805 0026

## OR5M11 Antibody [#G647]

**Catalog Number:** G647

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q96RB7

**Other Names:** olfactory receptor 5M11; olfactory receptor, family 5, subfamily M, member 11; OR11-199; OR5MB

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

OR5M11.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** OR5M11 antibody detects endogenous levels of total OR5M11 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:5000

**References:**

Fuchs T., Genomics 80:295-302(2002).

