



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR5P3 Antibody [#G652]

Catalog Number: G652

Concentration: 1mg/ml

Swiss-Prot No.: Q8WZ94

Other Names: JCG1; OL508; olfactory receptor 5P3;
olfactory receptor family 5 subfamily P member 3;
olfactory receptor OR11-94; olfactory receptor-like protein
JCG1

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from internal of human

OR5P3.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.

Specificity: OR5P3 antibody detects endogenous levels
of total OR5P3 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Gaudin J.-C., Chem. Senses 26:1157-1166(2001).

Suwa M., Submitted (JUL-2001) to the
EMBL/GenBank/DDBJ databases.

