



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

P2RY4 Antibody [#G708]

Catalog Number: G708

Concentration: 1mg/ml

Swiss-Prot No.: P51582

Other Names: P2Y purinoceptor 4; P2Y4; Uridine nucleotide receptor; UNR; P2P; P2RY4

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human P2RY4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: P2RY4 antibody detects endogenous levels of total P2RY4 protein.

Reactivity: Human, Mouse, Rat

Applications: WB: 1:500~1:1000 ELISA: 1:10000

References:

Communi D., J. Biol. Chem. 270:30849-30852(1995).

Nguyen T., J. Biol. Chem. 270:30845-30848(1995).

Stam N.J., FEBS Lett. 384:260-264(1996).

