



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

P2RY8 Antibody [#G711]

Catalog Number: G711

Concentration: 1mg/ml

Swiss-Prot No.: Q86VZ1

Other Names: G-protein coupled purinergic receptor P2Y8; P2Y purinoceptor 8; P2Y8; purinergic receptor P2Y, G-protein coupled, 8

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human P2RY8.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: P2RY8 antibody detects endogenous levels of total P2RY8 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Ross M.T., Nature 434:325-337(2005).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Adrian K., Biochim. Biophys. Acta 1492:127-138(2000).

