



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

P2RY11 Antibody [#G714]

Catalog Number: G714

Concentration: 1mg/ml

Swiss-Prot No.: Q96G91

Other Names: P2Y purinoceptor 11; P2Y11

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human P2RY11.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Specificity: P2RY11 antibody detects endogenous levels of total P2RY11 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Communi D.; Submitted (JAN-2001) to the EMBL/GenBank/DDBJ databases.

King M.M.; Submitted (OCT-2003) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team; Genome Res. 14:2121-2127(2004).

