



Web: [www.anbobio.com](http://www.anbobio.com)

[www.dijibio.com](http://www.dijibio.com)

Order: [order@dijibio.com](mailto:order@dijibio.com)

Tel: (86) 519 8805 0026

## P2RY13 Antibody [#G716]

**Catalog Number:** G716

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q9BPV8

**Other Names:** P2Y purinoceptor 13; P2Y13; G-protein coupled receptor 86; G-protein coupled receptor 94; GPR86, GPR94; FKSG77

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

P2RY13.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** P2RY13 antibody detects endogenous levels of total P2RY13 protein.

**Reactivity:** Human, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Wittenberger T., J. Mol. Biol. 307:799-813(2001).

Lee D.K., Gene 275:83-91(2001).

Communi D., J. Biol. Chem. 276:41479-41485(2001).

