



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

PTGER3 Antibody [#G721]

Catalog Number: G721

Concentration: 1mg/ml

Swiss-Prot No.: P43115

Other Names: Prostanoid EP3 receptor; PGE receptor, EP3 subtype; PGE2-R

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human PTGER3.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: PTGER3 antibody detects endogenous levels of total PTGER3 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:10000

References:

Adam M., FEBS Lett. 338:170-174(1994).

Schmid A., Eur. J. Biochem. 228:23-30(1995).

Yang J., Biochem. Biophys. Res. Commun. 198:999-1006(1994).

