



Web: [www.anbobio.com](http://www.anbobio.com)

[www.dijibio.com](http://www.dijibio.com)

Order: [order@dijibio.com](mailto:order@dijibio.com)

Tel: (86) 519 8805 0026

## PRLHR Antibody [#G728]

**Catalog Number:** G728

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P49683

**Other Names:** G protein-coupled receptor 10; GR3; PRLHR; prolactin-releasing peptide receptor; PrRP receptor; PrRPR

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human

PRLHR.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** PRLHR antibody detects endogenous levels of total PRLHR protein.

**Reactivity:** Human, Mouse, Rat

**Applications:** WB: 1:500~1:3000 ELISA: 1:20000

**References:**

Marchese A., Genomics 29:335-344(1995).

Hinuma S., Nature 393:272-276(1998).

Ozawa A., Mol. Endocrinol. 16:785-798(2002).

