



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR10AG1 Antibody [#G417]

Catalog Number: G417

Concentration: 1mg/ml

Swiss-Prot No.: Q8NH19

Other Names: O10AG; olfactory receptor 10AG1;
olfactory receptor OR11-160; olfactory receptor, family 10,
subfamily AG, member 1; OR11-160

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered
saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl,
0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against
synthesized peptide derived from internal of human
OR10AG1.

Purification: The antibody was affinity-purified from
rabbit antiserum by affinity-chromatography using
epitope-specific immunogen.


Specificity: OR10AG1 antibody detects endogenous
levels of total OR10AG1 protein.

Reactivity: Human


Applications: WB: 1:500~1:3000 IF: 1:100~1:500
ELISA: 1:40000

References:

Suwa M., Submitted (JUL-2001) to the
EMBL/GenBank/DDBJ databases.
Malnic B., Proc. Natl. Acad. Sci. U.S.A.
101:2584-2589(2004).



Western blot analysis of extracts from MCF-7 cells,
HUVEC cells and Jurkat cells, using OR10AG1 antibody.



Immunofluorescence analysis of MCF-7 cells, using
OR10AG1 antibody.