



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** +86 519 8805 0026

## OR10C1 Antibody [#G418]

**Catalog Number:** G418

**Concentration:** 1mg/ml

**Swiss-Prot No.:** Q96KK4

**Other Names:** Olfactory receptor 10C1; Hs6M1-17;  
OR10C1; OR10C2

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from internal of human OR10C1.

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** OR10C1 antibody detects endogenous levels of total OR10C1 protein.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:40000

**References:**

Mungall A.J., Nature 425:805-811(2003).

Fuchs T., Genomics 80:295-302(2002).

Menashe I., Nat. Genet. 34:143-144(2003).

