



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR2AK2 Antibody [#G434]

Catalog Number: G434

Concentration: 1mg/ml

Swiss-Prot No.: Q8NG84

Other Names: Olfactory receptor 2AK2; Olfactory receptor OR1-47; OR2AK2; OR2AK1P

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human OR2AK2.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR2AK2 antibody detects endogenous levels of total OR2AK2 protein.

Reactivity: Human (Identities = 100%, Positives = 100%)

Applications: IF: 1:100~1:500 ELISA: 1:5000

References:

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

Gregory S.G., Nature 441:315-321(2006).

Malnic B., Proc. Natl. Acad. Sci. U.S.A. 101:2584-2589(2004).

