



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: +86 519 8805 0026

OR2T2/2T35 Antibody [#G440]

Catalog Number: G440

Concentration: 1mg/ml

Swiss-Prot No.: Q8NGX2

Other Names: o2t35; olfactory receptor 2t35; olfactory receptor or1-66; olfactory receptor, family 2, subfamily t, member 35; or2t35; seven transmembrane helix receptor

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human OR2T2/2T35.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: OR2T2/2T35 antibody detects endogenous levels of total OR2T2/2T35 protein.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:1000

References:

Suwa M., Submitted (JUL-2001) to the EMBL/GenBank/DDBJ databases.

Malnic B., Proc. Natl. Acad. Sci. U.S.A. 101:2584-2589(2004).

