



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CATL1 (heavy chain, Cleaved-Thr288) Antibody [#L0184]

Catalog Number: L0184

Concentration: 1mg/ml

Swiss-Prot No.: P07711

Other Names: cathepsin L1; cathepsin L1 light chain; CATL; CATL1; CTSL; EC 3.4.22.15; FLJ31037; major excreted protein; MEP

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human CATL1.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CATL1 (heavy chain, Cleaved-Thr288) antibody detects endogenous levels of fragment of activated CATL1 resulting from cleavage adjacent to Thr288.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:10000

References:

Gal S., Biochem. J. 253:303-306(1988).

Joseph L.J., J. Clin. Invest. 81:1621-1629(1988).

Ebert L., Submitted (JUN-2004) to the EMBL/GenBank/DDBJ databases.

