



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CD97 β (Cleaved-Ser531) Antibody [#L0192]

Catalog Number: L0192

Concentration: 1mg/ml

Swiss-Prot No.: P48960

Other Names: CD97 antigen [Precursor], Leukocyte antigen CD97, CD97 antigen subunit alpha, CD97 antigen subunit beta, CD97

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human CD97 β .

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CD97 β (Cleaved-Ser531) antibody detects endogenous levels of fragment of activated CD97 β resulting from cleavage adjacent to Ser531.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:5000

References:

Hamann J., J. Immunol. 155:1942-1950(1995).

Hamann J., Genomics 32:144-147(1996).

Gray J.X., J. Immunol. 157:5438-5447(1996).

