



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

FA10 (light chain, Cleaved-Ala41) Antibody [#L0199]

Catalog Number: L0199

Concentration: 1mg/ml

Swiss-Prot No.: P00742

Other Names: Coagulation factor X [Precursor], EC 3.4.21.6, Stuart factor, Stuart-Prower factor, Factor X light chain, Factor X heavy chain, Activated factor Xa heavy chain, F10

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human FA10.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FA10 (light chain, Cleaved-Ala41) antibody detects endogenous levels of fragment of activated FA10 resulting from cleavage adjacent to Ala41.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Messier T.L., Gene 99:291-294(1991).

Leytus S.P., Biochemistry 25:5098-5102(1986).

Fung M.R., Proc. Natl. Acad. Sci. U.S.A. 82:3591-3595(1985).

