



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

FA12 (heavy chain, Cleaved-Arg372) Antibody [#L0203]

Catalog Number: L0203

Concentration: 1mg/ml

Swiss-Prot No.: P00748

Other Names: Coagulation factor XII [Precursor], EC 3.4.21.38, Hageman factor, HAF, Coagulation factor XIIIa heavy chain, Beta-factor XIIIa part 1, Beta-factor XIIIa part 2, Coagulation factor XIIIa light chain, F12

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human FA12.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: FA12 (heavy chain, Cleaved-Arg372) antibody detects endogenous levels of fragment of activated FA12 resulting from cleavage adjacent to Arg372.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Cool D.E., J. Biol. Chem. 262:13662-13673(1987).
Triptodi M., Nucleic Acids Res. 14:3146-3146(1986).
Cool D.E., J. Biol. Chem. 260:13666-13676(1985).

