



**Web:** [www.anbobio.com](http://www.anbobio.com)  
[www.dijibio.com](http://www.dijibio.com)  
**Order:** [order@dijibio.com](mailto:order@dijibio.com)  
**Tel:** (86) 519 8805 0026

## C1R (heavy chain, Cleaved-Arg463) Antibody [#L0229]

**Catalog Number:** L0229

**Concentration:** 1mg/ml

**Swiss-Prot No.:** P00736

**Other Names:** C1R; C1RA; Complement C1r component precursor

**Storage/Stability:** Store at -20°C/1 year

**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

**Immunogen:** The antiserum was produced against synthesized peptide derived from human C1R heavy chain.

**Purification:** The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Specificity:** C1R (heavy chain, Cleaved-Arg463) antibody detects endogenous levels of fragment of activated C1R heavy chain resulting from cleavage adjacent to Arg463.

**Reactivity:** Human

**Applications:** WB: 1:500~1:3000 ELISA: 1:10000

**References:**

Leytus S.P., Biochemistry 25:4855-4863(1986).

Journet A., Biochem. J. 240:783-787(1986).

Nakagawa M., Ann. Hum. Genet. 67:207-215(2003).

