



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CASP4 (p20, Cleaved-Gln81) Antibody [#L0156]

Catalog Number: L0156

Concentration: 1mg/ml

Swiss-Prot No.: P49662

Other Names: CASP-4; caspase 4, apoptosis-related cysteine peptidase; caspase-4; caspase-4 2; EC 3.4.22.57; EC 3.4.22.64; ICE(rel)-II; ICEREL-II; ICH-2; ICH-2 protease; ICH2; Mih1/TX; TX; TX protease

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human CASP4.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CASP4 (p20, Cleaved-Gln81) antibody detects endogenous levels of fragment of activated CASP4 resulting from cleavage adjacent to Gln81.

Reactivity: Human,

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Faucheu C., EMBO J. 14:1914-1922(1995).
Munday N.A., J. Biol. Chem. 270:15870-15876(1995).
Kamens J., J. Biol. Chem. 270:15250-15256(1995).

