



Web: www.anbobio.com

www.dijibio.com

Order: order@dijibio.com

Tel: (86) 519 8805 0026

Caspase 4/5 (p20, Cleaved-Asp270/Asp311) Antibody [#L0157]

Catalog Number: L0157

Concentration: 1mg/ml

Swiss-Prot No.: P49662/P51878

Other Names: DNA primase small subunit; DNA primase 49 kDa subunit; p49

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human Caspase 4/5.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: Caspase 4/5 (p20, Cleaved-Asp270/Asp311) antibody detects endogenous levels of fragment of activated Caspase 4/5 resulting from cleavage adjacent to Asp270/Asp311.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Faucheu C., EMBO J. 14:1914-1922(1995).

Munday N.A., J. Biol. Chem. 270:15870-15876(1995).

Kamens J., J. Biol. Chem. 270:15250-15256(1995).

