



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CASP5 (p10, Cleaved-Ser331) Antibody [#L0159]

Catalog Number: L0159

Concentration: 1mg/ml

Swiss-Prot No.: P51878

Other Names: Caspase-5 [Precursor], EC 3.4.22.58, CASP-5, ICH-3 protease, TY protease, ICE(rel)-III, Caspase-5 subunit p20, Caspase-5 subunit p10, CASP5, ICH3

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against

synthesized peptide derived from human CASP5.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CASP5 (p10, Cleaved-Ser331) antibody detects endogenous levels of fragment of activated CASP5 resulting from cleavage adjacent to Ser33.

Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Munday N.A., J. Biol. Chem. 270:15870-15876(1995).
Faucheu C., Eur. J. Biochem. 236:207-213(1996).

