



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CATD (light chain, Cleaved-Gln161) Antibody [#L0174]

Catalog Number: L0174

Concentration: 1mg/ml

Swiss-Prot No.: P07339

Other Names: CATD; Cathepsin D; Cathepsin D heavy chain; Cathepsin D light chain; ceroid-lipofuscinosis, neuronal 10; CLN10; CPSD; CTSD; lysosomal aspartyl peptidase; lysosomal aspartyl protease; MGC2311

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human CTSD.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CATD (light chain, Cleaved-Gln161) antibody detects endogenous levels of fragment of activated CTSD resulting from cleavage adjacent to Gln161.


Reactivity: Human

Applications: WB: 1:500~1:3000 ELISA: 1:20000

References:

Redecker B., DNA Cell Biol. 10:423-431(1991).

May F.E., Gene 134:277-282(1993).



For Research Use Only