



Web: www.anbobio.com
www.dijibio.com
Order: order@dijibio.com
Tel: (86) 519 8805 0026

CATD (heavy chain, Cleaved-Leu169) Antibody [#L0176]

Catalog Number: L0176

Concentration: 1mg/ml

Swiss-Prot No.: P07339

Other Names: Cathepsin D [Precursor], EC 3.4.23.5, Cathepsin D light chain, Cathepsin D heavy chain, CTSD, CPSD

Storage/Stability: Store at -20°C/1 year

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunogen: The antiserum was produced against synthesized peptide derived from human CATD.

Purification: The antibody was affinity-purified from

rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity: CATD (heavy chain, Cleaved-Leu169) antibody detects endogenous levels of fragment of activated CATD resulting from cleavage adjacent to Leu169.

Reactivity: Human


Applications: WB: 1:500~1:3000 ELISA: 1:40000

References:

Faust P.L., Proc. Natl. Acad. Sci. U.S.A. 82:4910-4914(1985).

Westley B.R., Nucleic Acids Res. 15:3773-3786(1987).

Redecker B., DNA Cell Biol. 10:423-431(1991).



For Research Use Only